

THE CUST

STRAUMANN® RN/WN synOcta® GOLD ABUTMENT

New

Straumann® RN synOcta® Gold Abutment

Straumann® WN synOcta® Gold Abutment

Intended use

- Cement-retained or screw-retained single crowns
- Telescopic constructions
- Cement-retained bridges on a mesostructure

Characteristics


Simple

- Reduced height for individual and esthetic results
- One-piece design for easy handling

Reliable

- Precise fit*

* Only Straumann® abutments on Straumann® implants ensure the most precise fit.





TECHNICAL INFORMATION

Components	RN and WN	
▪ Abutment	Ceramicor®*	
▪ Modelling aid	Burn-out plastic (POM)	
▪ Occlusal screw (Art. No. 048.356)	Titanium	
Specifications	RN	WN
▪ Abutment (height from implant shoulder)	4,3 mm	4,2 mm
▪ Modelling aid (height)	8,8 mm	8,8 mm
▪ Occlusal screw (length)	6,7 mm	6,7 mm
▪ Abutment with mounted modelling aid (length from abutment shoulder)	12,0 mm	11,9 mm
▪ Weight of gold abutment	ca. 0,61 g	ca. 0,78 g

* Ceramicor® is a registered trademark of Cendres & Métaux SA, Biel-Bienne (Switzerland).

ORDERING INFORMATION

Art. No.		Product
048.642		Straumann® RN synOcta® Gold Abutment*
048.644		Straumann® WN synOcta® Gold Abutment*

* includes basal screw (also separately available as Art. No. 048.356).



STRAUMANN
GUARANTEE

International Headquarters

Institut Straumann AG
 Peter Merian-Weg 12
 CH-4002 Basel, Switzerland
 Phone +41 (0)61 965 11 11
 Fax +41 (0)61 965 11 01